

Enduras WIRELESS 10 x 12



Enduras Wireless 10 x 12, MinXray's most complete, portable and rugged digital radiography system is now completely wireless. This compact, battery-operated digital system delivers superior image performance while easily transported to the patient. Features of the **Enduras Wireless 10 x 12** include:

- Cassette-sized DR panel—will fit conventional bucky tray
- Auto-sensing panel compatible with any x-ray generator
- Exceptional image quality with CsI Scintillator

MINXRAY[®] INC

Where you care, we're there.

Enduras WIRELESS 10 x 12

CsI FLAT PANEL DETECTOR SPECIFICATIONS

SCINTILLATOR	Direct Deposit Cesium Iodide (CsI)
ACTIVE AREA	10 x 12 in. (24.4 x 30.7 cm)
PIXEL ARRAY	2048 x 2560
PIXEL PITCH	120 µm
GRayscale	16 bits
IMAGE ACQUISITION TIME	4 - 5 seconds
EXPOSURE CONTROL	Full Field Automatic Exposure Detection (F ² AED) or manual
DIMENSION	11 x 13 x 0.6 in. (28.2 x 33.3 x 1.5 cm)
WEIGHT	5.1 lbs (2.3 kg) <i>including battery</i>
HOUSING MATERIAL	Carbon Fiber Front - High Strength Aluminum Alloy Back
OPERATING TEMPERATURE RANGE	40 - 95° F (5 - 35° C)
SUPPLY	100 - 250V AC
FREQUENCY	50/60 Hz
MANUFACTURER'S WARRANTY	5 years

X-RAY UNIT

	AVAILABLE kVDC	AVAILABLE MA	UNIT WEIGHT
TR8020 ULTRA LIGHT	50-80	20-16	14.5 lb (6.3 kg)
TR90 TOUGH RAY	40-90	30-9	14.1 lb (6.4 kg)
TR90B BATTERY POWERED	40-90	20-10	15.0 lb (6.8 kg)
HF100+	40-100	40-20	24.2 lb (11 kg)



Shown with optional TR90B

Ordering Information

Enduras Wireless 10 x 12 Portable DR System

- 1 ea. CsI Flat panel detector
- 2 ea. Lithium ion rechargeable batteries
- 1 ea. Laptop or tablet
- 1 ea. MinXView Imaging Software
- 1 ea. Exposure switch & cord, power cord & carrying case

Suggested Optional Portable MinXray Units:

- 1 ea. **TR8020**
- 1 ea. **TR90**
- 1 ea. **TR90B**
- 1 ea. **HF100+**

For more information about MinXray products, contact your local dealer or MinXray at:

MinXray, Inc.
3611 Commercial Avenue
Northbrook, Illinois 60062 USA
1-847-564-0323

Toll Free: 1-800-221-2245 (U.S. & CA.)
E-mail: info@minxray.com
Web: www.minxray.com

MINXRAY[®] INC

Where you care, we're there.